



HUMAN HEALTH

ENVIRONMENTAL HEALTH


ACCURATE.
BRIGHTER.
BETTER.

Distinctive Luminescence
Instruments and Reagents



PerkinElmer[®]
For the Better

TABLE OF CONTENTS



LUMINESCENCE ASSAY SYSTEMS AND INSTRUMENTS

REPORTER GENE ASSAYS

| | |
|--|---|
| neolite™ | 2 |
| britelite™ plus and steadylite plus™ | 3 |

ATP QUANTIFICATION

| | |
|---|---|
| ATPlite® cytotoxicity and cell proliferation assays | 5 |
|---|---|

CALCIUM FLUX ASSAY SYSTEMS

| | |
|--|---|
| AequoScreen® and PhotoScreen® assays | 7 |
|--|---|

CALCIUM DETECTION

| | |
|---|---|
| AequoZen® ready-to-use frozen cells | 8 |
|---|---|

CELL LINES

| | |
|----------------------------|---|
| Ordering Information | 9 |
|----------------------------|---|

DETECTION INSTRUMENTATION AND AUTOMATION

| | |
|--|----|
| EnSpire™ Multilabel Plate Reader with Ultra-sensitive Luminescence | 12 |
| VICTOR™ Plate Readers | 12 |
| EnVision® Multilabel Readers | 12 |
| MicroBeta ² ® and TopCount® Plate Counters | 13 |
| ViewLux™ Imagers | 13 |
| JANUS® Automated Workstation | 13 |

INSTRUMENT QUICK REFERENCE GUIDE

| | |
|---|----|
| Microplates for Luminescence Assays | 16 |
|---|----|



STREAMLINED DETECTION. SIMPLY BRILLIANT.

Luminescence detection platforms provide simple and reproducible results over a wide dynamic range offering great value for many research and discovery applications. Choose from a variety of assay systems and instruments to optimally streamline your detection needs.

KEY BENEFITS

Luminescence detection:

- **Cost-effective**
- **Easy set-up** for high throughput
- **Rapid** – miniaturizable and automatable; some systems allow results in 15-25 minutes
- **No separation** steps (homogeneous assays)
- **Sensitive** – better than fluorescence technologies
- **Less interference**

LUMINESCENCE ASSAY SYSTEMS AND INSTRUMENTS

| | AequoScreen AequoZen PhotoScreen | britelite plus | steadylite plus | neolite | ATPlite 1step | ATPlite |
|-----------------------------|--|--|--|--|--|--|
| Applications | <ul style="list-style-type: none"> • Calcium Flux • GPCR • Ion Channels | <ul style="list-style-type: none"> • Reporter Gene Assay • GPCR | <ul style="list-style-type: none"> • Reporter Gene Assay • GPCR | <ul style="list-style-type: none"> • Reporter Gene Assay • GPCR | <ul style="list-style-type: none"> • Cytotoxicity • Cell Proliferation Assay | <ul style="list-style-type: none"> • Cytotoxicity • Cell Proliferation Assay |
| Half-life (hours) | Flash | 0.5 | 4 to 5 | 2.5 | 0.5 | 5 |
| Sensitivity | High | Very High | Moderate | High | As few as 5 cells/well | As few as 5 cells/well |
| Plate Format (wells) | 96 384 1536 | 96 384 1536 | 96 384 1536 | 96 384 1536 | 96 384 1536 | 96 384 1536 |
| Instrumentation | — MicroBeta ² LumiJET EnVision VICTOR X VICTOR X Light — | TopCount MicroBeta ² EnSpire EnVision VICTOR X VICTOR X Light ViewLux | TopCount MicroBeta ² EnSpire EnVision VICTOR X VICTOR X Light ViewLux | TopCount MicroBeta ² EnSpire EnVision VICTOR X VICTOR X Light ViewLux | TopCount MicroBeta ² EnSpire EnVision VICTOR X VICTOR X Light ViewLux | TopCount MicroBeta ² EnSpire EnVision VICTOR X VICTOR X Light ViewLux |



INTENSITY. FLEXIBILITY. JUST RIGHT FOR EVERY ASSAY.

Whether your assay requires the sensitivity provided by high signal intensity, or the flexibility of an extended signal half-life, PerkinElmer provides a solution. Choose from one of three luciferase detection systems: neolite, britelite plus, or steadylite plus depending on the specific requirements of your experiment.


KEY BENEFITS

- **No more washing** – just add your lites reagent of choice, mix, incubate, and read
- **Hassle-free storage** – all lites reagents are designed for extended storage at 2-8 °C (no -20 °C storage required)
- **Odor-friendly**, eliminates hood work – lites products do not contain DTT and do not have the hazards and repugnant odor typical of conventional luciferase detection systems
- **Excellent Z' values** – superior S/B ratios allow for excellent Z' values in antagonist and agonist studies

neolite — Bright and Stable

neolite offers the ultimate compromise between strong signal intensity and extended signal half-life. It provides a minimum half-life of 2.5 hours while maintaining excellent sensitivity. neolite also delivers superior reproducibility, as it has less sensitive

mixing conditions than typical luminescent detection systems. It is ideal for medium to high throughput applications which employ batch processing methods but still require exceptional sensitivity.



steadylite plus — Long-lived Glow


steadylite plus was developed to provide an extended signal half-life of up to 5 hours. This extended half-life makes steadylite plus the reagent of choice for high throughput screening where you require the same steady signal first plate in to last plate out.

britelite plus — Intensively Sensitive

britelite plus possesses an extremely strong signal intensity. With a signal half-life of 30 minutes, it is ideal for assays requiring the utmost sensitivity and for any application employing continuous processing methodology. See Applications Table on page 4.




EnSpire multilabel plate reader with ultra-sensitive luminescence requires fewer cells to get the same signal as standard luminescence saving on expensive transfection reagent and substrate costs.



A dilution series of luciferase enzyme where luciferase, britelite plus and neolite have a sub-femtomole detection limit for luciferase compared to steadylite plus.

Luminescence for 10⁻¹⁰ g luciferase per well



Time course of the luminescent signal of britelite plus, neolite and steadylite plus where britelite plus shows an intense signal with rapid decrease in signal and neolite and steadylite plus both show strong signal intensity over an extended period of time.





The PerkinElmer suite of luciferase detection systems brings the simplicity of reporter gene assays to a whole new level. All lites

products offer superior performance, while still providing convenient one-step processing and ease of storage and handling.

| | neolite | steadylite plus | britelite plus |
|-------------------|--|--|--|
| Sensitivity | High | Moderate | Very High |
| Half-life (hours) | 2.5 | 5 | 0.5 |
| Plate Format | 96, 384, 1536 | 96, 384, 1536 | 96, 384, 1536 |
| Applications | <ul style="list-style-type: none"> • Low transfection efficiencies • Primary cell transfection • Increased assay windows • Extended read times • Batch processing | <ul style="list-style-type: none"> • Strong steady signal • Extended read time • High throughput screening • Extended batch processing | <ul style="list-style-type: none"> • Superior sensitivity • Low transfection efficiencies • Stem cell transfection • Rapid read time – continuous processing |

| Part Number | | | Volume (mL) | Data Points* in 96-, 384- and 1536-well microplates |
|-------------|-----------------|----------------|-------------|--|
| neolite | steadylite plus | britelite plus | | |
| 6016716 | 6016756 | 6016766 | 10 | 100, 400, 3300 |
| 6016711 | 6016751 | 6016761 | 100 | 1000, 4000, 33000 |
| 6016767 | 6016757 | 6016767 | 500 | 5000, 20000, 165000 |
| 6016769 | 6016759 | 6016769 | 1000 | 10000, 40000, 330000 |

* The recommended assay volumes are 100 μ L for 96-well microplates, 25 μ L for 384-well microplates and 3 μ L for 1536-well microplates.



Agonist and antagonist dose-response curves comparing britelite plus, neolite and steadylite plus.



LIGHT UP THE ATP LEVELS AT YOUR PACE.

Fast, reliable detection of ATP for cell quantification

ATPlite 1step and ATPlite are patented¹ innovative technologies that measure cell proliferation and cytotoxicity in mammalian cells based on the

detection of ATP using firefly luciferase. Light production caused by the reaction of ATP with added luciferase and D-luciferin is proportional to the ATP concentration.

ATP is a marker for cell viability because it is present in all metabolically active cells. Because ATP concentration declines rapidly when cells undergo necrosis or apoptosis, monitoring ATP is a good indicator of cytotoxic, cytostatic, and proliferation effects.

KEY BENEFITS

ATPlite 1step – Highest Sensitivity

- True “mix and measure” continuous processing
- Ultra high sensitivity, 3 times more light output than ATPlite
- Short equilibrium time for faster processing
- Signal half-life of 30 minutes

ATPlite – Convenient, Extended Signal

- 2-step format, non-separation assay for time course studies; lysate can be stored frozen and tested later
- Quantitate from <5 cells/well
- Prolonged equilibrium time for batch processing
- Signal half-life of 5 hours

| Description | Part Number | Volume (mL) | Data Points* in 96-, 384- and 1536-well microplates |
|---------------|-------------|-------------|--|
| ATPlite 1step | 6016736 | 10 | 100, 400, 3300 |
| | 6016731 | 100 | 1000, 4000, 33000 |
| | 6016739 | 1000 | 10000, 40000, 330000 |
| ATPlite | 6016943 | 30 | 300, 1200, 9900 |
| | 6016941 | 100 | 1000, 4000, 33000 |
| | 6016947 | 500 | 5000, 20000, 16500 |
| | 6016949 | 1000 | 10000, 40000, 330000 |

* The recommended assay volumes are 100 µL for 96-well microplates, 25 µL for 384-well microplates and 3 µL for 1536-well microplates.

1. U.S. Pat. 6503723; EP Pat. 117825 (CH, DK, GB, DE 69924127);
NL Pat. 1010224; Canada Patent Appl. No. 2345721; Australian Patent Appl. No. 754602; other patents pending.



ROBUST ASSAYS FOR CALCIUM FLUX IN A FLASH.

AequoScreen Aequorin and PhotoScreen Photina photoproteins are derived from the jellyfish *Aequoria victoria* and *Obelia longissima* respectively. Upon calcium binding, the photoprotein oxidizes coelenterazine into coelenteramide with the production of CO₂ and emission of light. Because both can reduce costs and increase productivity, they are ideal replacement technologies for standard fluorescent dyes.


KEY BENEFITS:

Photoprotein technology:


- **Fully automatable** – can be miniaturized to 1536 format
- **Reduces reagent and consumable costs**, including the number of cells required for an assay
- **Allows accurate detection** of agonists, antagonists, partial agonists, and allosteric modulators through high signal-to-noise background ratio
- **Delivers robust assay** (high Z', low CV values)
- **Offers powerful tools** to convert mammalian cell lines already expressing a GPCR to a luminescence assay format
- **Reduces false positives** from auto-fluorescent compounds in screening libraries
- **Increases throughput** using option to work with cells in suspension
- **Achieves efficient, high throughput** – over 200,000 data points/day (depending on assay instrument and microplate format)

The ideal calcium flux assay for GPCR analysis — AequoScreen

AequoScreen is a generic GPCR technology, which can be used with dedicated flash luminescence instrumentation (MicroBeta² LumiJET, VICTOR, EnVision, FLIPR[®]TETRA, FLIPR[®]TETRA (aequorin option) FDSS) to measure calcium from G_q-coupled receptors or G_s, G_i-coupled receptors through a promiscuous or chimeric G-protein.



Histamine H₁ AequoScreen cell line agonist and antagonist suspension assay dose response curves demonstrate the quality of data produced by the MicroBeta² LumiJET. Histamine H₁ AequoScreen cells (5000 cells) were dispensed together with buffer or antagonist (trans-triprolidine, final concentration 10 x IC₅₀) into an OptiPlate-384 in 36 replicates. The agonist (histamine, final concentration 10 x EC₅₀), was dispensed and resulting luminescence was recorded using the MicroBeta² LumiJET (total volume, 70 μl).



Histamine dose response curve, 1536-well suspension assay, 1500 cells/well (mean +/- SD from 32 wells) run on a LumiLux[®], demonstrating that the aequorin technology can be miniaturized.

The AequoScreen platform includes:


- **Parental cell lines** expressing mitochondrially targeted apoaequorin with or without the promiscuous G-protein $G\alpha_{16}$
- **Double transfected cell lines** expressing mitochondrially targeted apoaequorin and GPCR with or without the promiscuous G-protein $G\alpha_{16}$
- **AequoZen** frozen ready-to-use cells: growth arrested double transfected cell lines available frozen for immediate use (see next page)
- **Aequorin plasmids** expressing mitochondrially targeted apoaequorin protein under the control of the CMV promoter, with or without the promiscuous G-protein $G\alpha_{16}$
- **Fully validated Starter Kits** for AequoScreen cell lines and AequoZen ready-to-use cells enabling you to become familiar with the AequoScreen technology

PhotoScreen — a novel assay for calcium coupled ion channel analysis


Photina is a chimeric photoprotein, generated from obelin and clytin, which can be used to measure calcium from G_T , G_q , and G_s - coupled receptors as well as calcium-coupled ion channels and transporters. Photina light emission can also be measured using flash luminescence instrumentation.

The PhotoScreen platform includes:


- **Double transfected cell lines** expressing mitochondrially targeted apophotoprotein and a recombinant calcium coupled ion channel/transporter or a recombinant GPCR with or without the $G_{q/3}$
- **Fully Validated Cell Line Starter Kit**



IB-MECA agonist 384-well adherent assay, cell density study. (Mean +/- SD 8 wells) run on a FLIPR[®]TETRA.



The Muscarinic M_5 AequoScreen cell line agonist (Acetylcholine) and antagonist (N-Me-Scopolamine) dose response curves run on a MicroBeta² LumiJET. 384-well format, 6,000 cells/well.



384-well suspension assay, cell density study, (Mean +/- SD 8 wells) run on a LumiLux.



FULLY VALIDATED FROZEN CELLS OFF THE SHELF OR CUSTOM-MADE

AequoZen frozen cells, pre-validated for aequorin calcium detection, enable rapid and cost-effective GPCR profiling, pharmacology or screening applications.

KEY BENEFITS

- **Available off-the-shelf** or customized to your needs
- **Simple** step-by-step protocols for cell culture functional testing
- **Affordable** access to small aliquots of cells (10 x 10⁶ cells)
- **Highest quality**, growth-arrested gamma irradiated cells to minimize biological variation
- **Versatile** up-scaling to suit specific needs

Ready-to-use cells take the stress out of functional testing

We remove the lengthy cell culture process from your functional testing; we culture the cells, freeze them with an optimized protocol, and gamma irradiate them to stop replication. These cells are then stored in liquid nitrogen or at -80 °C. Simply thaw them and use in your functional, cellular GPCR test.

Confident validation



We validate defined families of GPCRs that let you perform selectivity studies rapidly and cost-effectively as needed. Validated batches give results meeting our QC criteria, both on agonist EC₅₀ and signal window, so you can obtain the same reliable results with every purchase, every time.

Convenient flexibility

Frozen cells are available off-the-shelf in aliquots of 10 million cells per vial only for catalog products. Other sizes are offered for custom service only. We also scale-up to provide you with validated cells for full screening campaigns. Whether it is for screening, lead optimization or profiling, you can now perform cellular GPCR tests immediately by using batches of frozen cells on multiple receptors at a time. Users have the option to work in both suspension and adherent mode.

Quality you can trust

Our frozen cell portfolio undergoes a rigorous QC procedure; this procedure does not leave any trace of hazardous chemicals in the cell media and has been titrated to preserve functional response.



All batches are tested for quality of assay window and consistency of EC₅₀ values. The calcium aequorin assay (left) illustrates the correlation between EC₅₀ values of the reference agonist obtained on fresh versus frozen cells with an aequorin readout on 15 receptors. The AT₁ pharmacology comparison values (right) are similar between fresh AequoScreen AT₁ cell line kept in continuous culture and equivalent frozen cells.

CELL LINE AND FROZEN CELL PRODUCT LISTING



ORDERING INFORMATION

If you do not see the target you are interested in, please inquire as we are continuously expanding our portfolio and also provide custom cell line development.

| STARTER KITS | | | |
|----------------------------|--|-----------------------|------------|
| Product | Receptor | Type | Cat. No. |
| AequoScreen starter kit | H ₁ and M ₅ | AequoScreen Cell Line | ES-001-A |
| AequoZen starter kit | H ₁ and M ₅ | AequoZen Frozen Cells | ES-001-AF |
| PhotoScreen starter kit | A ₃ | PhotoScreen Cell Line | AX-001-PCF |
| GPCR CELL LINES | | | |
| Receptor type | Subtype | Product Type | Cat. No. |
| 5-Hydroxytryptamine | 5-HT _{1A} | AequoScreen Cell Line | ES-310-A |
| | 5-HT _{1A} | AequoZen Frozen Cells | ES-310-AF |
| | 5-HT _{2A} | AequoScreen Cell Line | ES-313-A |
| | 5-HT _{2A} | AequoZen Frozen Cells | ES-313-AF |
| | 5-HT _{2B} | AequoScreen Cell Line | ES-314-A |
| | 5-HT _{2C} edited isoform | AequoScreen Cell Line | ES-315-A |
| | 5-HT _{2C} edited isoform | AequoZen Frozen Cells | ES-315-AF |
| | 5-HT _{2C} non edited isoform | AequoScreen Cell Line | ES-318-A |
| | 5-HT _{2C} non edited isoform | AequoZen Frozen Cells | ES-318-AF |
| | 5-HT _{3A} | AequoScreen Cell Line | ES-402-A |
| | 5-HT _{3A} | AequoZen Frozen Cells | ES-402-AF |
| | 5-HT _{5A} | AequoScreen Cell Line | ES-401-A |
| | 5-HT _{5A} | AequoZen Frozen Cells | ES-401-AF |
| | 5-HT ₆ | AequoScreen Cell Line | ES-316-A |
| 5-HT ₇ | AequoZen Frozen Cells | ES-316-AF | |
| Acetylcholine (Muscarinic) | M ₁ | AequoScreen Cell Line | ES-210-A |
| | M ₁ | AequoZen Frozen Cells | ES-210-AF |
| | M ₂ | AequoScreen Cell Line | ES-211-A |
| | M ₂ | AequoZen Frozen Cells | ES-211-AF |
| | M ₃ | AequoScreen Cell Line | ES-212-A |
| | M ₃ | AequoZen Frozen Cells | ES-212-AF |
| | M ₄ | AequoScreen Cell Line | ES-213-A |
| | M ₄ | AequoZen Frozen Cells | ES-213-AF |
| | M ₅ | AequoScreen Cell Line | ES-214-A |
| Adenosine | A ₁ | AequoScreen Cell Line | ES-010-A |
| | A ₁ | AequoZen Frozen Cells | ES-010-AF |

| Receptor type | Subtype | Product Type | Cat. No. |
|-------------------|-----------------------|-----------------------|------------|
| | A _{2A} | AequoScreen Cell Line | ES-011-A |
| | A _{2A} | AequoZen Frozen Cells | ES-011-AF |
| | A _{2B} | AequoScreen Cell Line | ES-013-A |
| | A ₃ | AequoScreen Cell Line | ES-012-A |
| | A ₃ | AequoZen Frozen Cells | ES-012-AF |
| | A ₃ | PhotoScreen Cell Line | AX-001-PCF |
| Adrenoceptors | alpha _{1A} | AequoScreen Cell Line | ES-036-A |
| | alpha _{1A} | AequoZen Frozen Cells | ES-036-AF |
| | alpha _{1B} | AequoScreen Cell Line | ES-037-A |
| | alpha _{1D} | AequoScreen Cell Line | ES-038-A |
| | alpha _{2A} | AequoZen Frozen Cells | ES-030-AF |
| | alpha _{2A} | AequoScreen Cell Line | ES-030-A |
| | alpha _{2B} | AequoScreen Cell Line | ES-031-A |
| | alpha _{2B} | AequoZen Frozen Cells | ES-031-AF |
| | alpha _{2C} | AequoScreen Cell Line | ES-032-A |
| | beta ₁ | AequoScreen Cell Line | ES-033-A |
| beta ₁ | AequoZen Frozen Cells | ES-033-AF | |
| Anaphylatoxin | C3a | AequoScreen Cell Line | ES-730-A |
| | C5a | AequoScreen Cell Line | ES-731-A |
| | C5a | AequoZen Frozen Cells | ES-731-AF |
| Angiotensin | AT ₁ | AequoScreen Cell Line | ES-072-A |
| | AT ₁ | AequoZen Frozen Cells | ES-072-AF |
| Apelin | APJ | AequoScreen Cell Line | ES-460-A |
| Beta Alanine | MRGPRD, TGR7 | AequoScreen Cell Line | ES-741-A |
| Bombesin | BB ₁ | AequoScreen Cell Line | ES-581-A |
| | BB ₂ | AequoScreen Cell Line | ES-582-A |
| | BB ₃ | AequoScreen Cell Line | ES-580-A |
| Bradykinin | B ₁ | AequoScreen Cell Line | ES-091-A |
| | B ₁ (rat) | AequoScreen Cell Line | ES-092-A |
| | B ₂ | AequoScreen Cell Line | ES-090-A |
| Calcitonin | AM ₂ | AequoScreen Cell Line | ES-430-A |
| | CGRP ₁ | AequoScreen Cell Line | ES-420-A |
| | CGRP ₁ | AequoZen Frozen Cells | ES-420-AF |
| Cannabinoid | CB ₁ | AequoScreen Cell Line | ES-110-A |
| | CB ₁ | AequoZen Frozen Cells | ES-110-AF |
| Chemerin | GPR1 | AequoScreen Cell Line | ES-661-A |

| Receptor type | Subtype | Product Type | Cat. No. |
|--------------------------------|-----------------------|-----------------------|-----------|
| Chemerin | GPR1 | AequoZen Frozen Cells | ES-661-AF |
| Chemokine | CCR10 | AequoScreen Cell Line | ES-143-A |
| | CCR10 | AequoZen Frozen Cells | ES-143-AF |
| | CCR2b | AequoScreen Cell Line | ES-133-A |
| | CCR2b | AequoZen Frozen Cells | ES-133-AF |
| | CCR3 | AequoScreen Cell Line | ES-138-A |
| | CCR3 | AequoZen Frozen Cells | ES-138-AF |
| | CCR6 | AequoScreen Cell Line | ES-139-A |
| | CCR6 | AequoZen Frozen Cells | ES-139-AF |
| | CCR7 | AequoScreen Cell Line | ES-140-A |
| | CCR7 | AequoZen Frozen Cells | ES-140-AF |
| | CCR8 | AequoScreen Cell Line | ES-136-A |
| | CCR8 | AequoZen Frozen Cells | ES-136-AF |
| | CCR9a | AequoScreen Cell Line | ES-146-A |
| | CCR9a | AequoZen Frozen Cells | ES-146-AF |
| | CX ₃ CR1 | AequoScreen Cell Line | ES-137-A |
| | CX ₃ CR1 | AequoZen Frozen Cells | ES-137-AF |
| | CXCR2 | AequoScreen Cell Line | ES-145-A |
| | CXCR2 | AequoZen Frozen Cells | ES-145-AF |
| | CXCR3 | AequoScreen Cell Line | ES-142-A |
| | CXCR3 | AequoZen Frozen Cells | ES-142-AF |
| | CXCR6 | AequoScreen Cell Line | ES-720-A |
| | CXCR6 | AequoZen Frozen Cells | ES-720-AF |
| | XCR1 | AequoScreen Cell Line | ES-148-A |
| XCR1 | AequoZen Frozen Cells | ES-148-AF | |
| CCR1 | AequoScreen Cell Line | ES-132-A | |
| Cholecystokinin | CCK ₁ | AequoScreen Cell Line | ES-530-A |
| | CCK ₁ | AequoZen Frozen Cells | ES-530-AF |
| | CCK ₂ | AequoScreen Cell Line | ES-531-A |
| | CCK ₂ | AequoZen Frozen Cells | ES-531-AF |
| Corticotropin-Releasing Factor | CRF ₁ | AequoScreen Cell Line | ES-152-A |
| Dopamine | D _{2L} | AequoScreen Cell Line | ES-171-A |
| | D _{2L} | AequoZen Frozen Cells | ES-171-AF |
| Endothelin | ET _A | AequoScreen Cell Line | ES-320-A |
| | ET _A | AequoZen Frozen Cells | ES-320-AF |
| | ET _B | AequoScreen Cell Line | ES-321-A |
| | ET _B | AequoZen Frozen Cells | ES-321-AF |
| Formyl Peptide | FPRL1 | AequoScreen Cell Line | ES-610-A |
| Free Fatty Acid | FFA1 (GPR40) | AequoScreen Cell Line | ES-652-A |
| | FFA1 (GPR40) | AequoZen Frozen Cells | ES-652-AF |
| | GPR120 | AequoScreen Cell Line | ES-800-A |
| | GPR120 | AequoZen Frozen Cells | ES-800-AF |

| Receptor type | Subtype | Product Type | Cat. No. |
|--------------------------------------|--------------------------------------|-----------------------|------------|
| GABA _B | GABA _{B1aB2} | AequoScreen Cell Line | ES-500-A |
| Galanin | GAL ₁ | AequoScreen Cell Line | ES-510-A |
| | GAL ₂ | AequoScreen Cell Line | ES-511-A |
| | GAL ₂ | AequoZen Frozen Cells | ES-511-AF |
| | GAL ₂ | PhotoScreen Cell Line | AX-008-PCF |
| Ghrelin | ghrelin | AequoScreen Cell Line | ES-410-A |
| | ghrelin | AequoZen Frozen Cells | ES-410-AF |
| Glucagon | glucagon | AequoScreen Cell Line | ES-710-A |
| | Secretin | AequoScreen Cell Line | ES-712-A |
| | Secretin | AequoZen Frozen Cells | ES-712-AF |
| Glycoprotein Hormone | TSH | AequoScreen Cell Line | ES-790-A |
| Gonadotropin-Releasing Hormone | GnRH | AequoScreen Cell Line | ES-600-A |
| Histamine | H ₁ | AequoScreen Cell Line | ES-390-A |
| | H ₁ | AequoZen Frozen Cells | ES-390-AF |
| | H ₂ | AequoScreen Cell Line | ES-391-A |
| | H ₂ | AequoZen Frozen Cells | ES-391-AF |
| | H ₃ | AequoScreen Cell Line | ES-392-A |
| | H ₃ | AequoZen Frozen Cells | ES-392-AF |
| | H ₄ | AequoScreen Cell Line | ES-393-A |
| | H ₄ | AequoZen Frozen Cells | ES-393-AF |
| KISS1 (Metastin) | KISS1 (GPR54) | AequoScreen Cell Line | ES-630-A |
| Leukotriene | BLT ₁ (LTB4R1) | AequoScreen Cell Line | ES-340-A |
| | CysLT1 | AequoScreen Cell Line | ES-470-A |
| | OXE (HM74-like) | AequoScreen Cell Line | ES-640-A |
| Lysophospholipid | S ₁ P ₂ (EDG5) | AequoScreen Cell Line | ES-594-A |
| | S ₁ P ₂ (EDG5) | AequoZen Frozen Cells | ES-594-AF |
| | S ₁ P ₄ (EDG6) | AequoScreen Cell Line | ES-592-A |
| | S ₁ P ₄ (EDG6) | AequoZen Frozen Cells | ES-592-AF |
| | S ₁ P ₅ (EDG8) | AequoScreen Cell Line | ES-593-A |
| S ₁ P ₅ (EDG8) | AequoZen Frozen Cells | ES-593-AF | |
| Mas related | MrgX1 | AequoScreen Cell Line | ES-740-A |
| | MrgX2 | AequoScreen Cell Line | ES-742-A |
| Melanin-Concentrating Hormone | MCH ₁ | AequoScreen Cell Line | ES-370-A |
| | MCH ₁ | AequoZen Frozen Cells | ES-370-AF |
| | MCH ₂ | AequoScreen Cell Line | ES-371-A |
| | MCH ₂ | AequoZen Frozen Cells | ES-371-AF |
| Melanocortin | MC ₄ | AequoScreen Cell Line | ES-191-A |
| | MC ₄ | AequoZen Frozen Cells | ES-191-AF |
| Melatonin | MT ₁ | AequoScreen Cell Line | ES-620-A |
| | MT ₂ | AequoScreen Cell Line | ES-621-A |
| | MT ₂ | AequoZen Frozen Cells | ES-621-AF |
| Motilin | motilin | AequoScreen Cell Line | ES-380-A |

| Receptor type | Subtype | Product Type | Cat. No. |
|-----------------------------|--------------------------------------|-----------------------|-----------|
| Neurokinin | NK ₂ | AequoScreen Cell Line | ES-251-A |
| | NK ₃ | AequoScreen Cell Line | ES-252-A |
| Neuromedin U | NMU1 | AequoScreen Cell Line | ES-450-A |
| | NMU2 | AequoScreen Cell Line | ES-451-A |
| Neuropeptide FF | NPF2 | AequoScreen Cell Line | ES-490-A |
| Neuropeptide S | NPS | AequoScreen Cell Line | ES-770-A |
| Neuropeptide Y | Y ₁ | AequoScreen Cell Line | ES-351-A |
| | Y ₂ | AequoScreen Cell Line | ES-352-A |
| | Y ₂ | AequoZen Frozen Cells | ES-352-AF |
| Neurotensin | NTS ₁ | AequoScreen Cell Line | ES-690-A |
| | NTS ₁ | AequoZen Frozen Cells | ES-690-AF |
| Nicotinic | GPR109A (HM74A) | AequoScreen Cell Line | ES-760-A |
| Opioid | kappa | AequoScreen Cell Line | ES-541-A |
| | mu | AequoScreen Cell Line | ES-542-A |
| | mu | AequoZen Frozen Cells | ES-542-AF |
| | NOP (ORL1) | AequoScreen Cell Line | ES-230-A |
| Orexin | OX ₁ | AequoScreen Cell Line | ES-330-A |
| | OX ₁ | AequoZen Frozen Cells | ES-330-AF |
| | OX ₂ | AequoScreen Cell Line | ES-331-A |
| | OX ₂ | AequoZen Frozen Cells | ES-331-AF |
| Oxoglutarate Receptor | OXGR1 (GPR99, Citric acid cycle) | AequoScreen Cell Line | ES-743-A |
| | OXGR1 (GPR99, Citric acid cycle) | AequoZen Frozen Cells | ES-743-AF |
| Peptide P518 | QRFP (OX ₂ -like, GPR103) | AequoScreen Cell Line | ES-670-A |
| Prokinectin | PKR ₁ | AequoScreen Cell Line | ES-750-A |
| | PKR ₁ | AequoZen Frozen Cells | ES-750-AF |
| | PKR ₂ | AequoScreen Cell Line | ES-751-A |
| Prolactin Releasing Peptide | PRRP | AequoScreen Cell Line | ES-480-A |
| Prostanoid | DP ₁ | AequoScreen Cell Line | ES-560-A |
| | DP ₂ (CRTH2) | AequoScreen Cell Line | ES-561-A |
| Protease-activated | PAR2 | AequoScreen Cell Line | ES-780-A |
| | PAR4 | AequoScreen Cell Line | ES-782-A |
| | PAR4 | AequoZen Frozen Cells | ES-782-AF |
| Purinergic | P2Y ₁₁ | AequoScreen Cell Line | ES-293-A |
| | P2Y ₁₁ | AequoZen Frozen Cells | ES-293-AF |
| Relaxin | RXFP3 | AequoScreen Cell Line | ES-656-A |
| | RXFP4 | AequoScreen Cell Line | ES-659-A |

| Receptor type | Subtype | Product Type | Cat. No. |
|-------------------------------|--------------------------|-----------------------|-----------------------|
| Somatostatin | sst _{2a} | AequoScreen Cell Line | ES-521-A |
| | sst ₄ | AequoZen Frozen Cells | ES-524-AF |
| | sst ₄ | AequoScreen Cell Line | ES-524-A |
| | sst ₅ | AequoScreen Cell Line | ES-522-A |
| Succinate | SUCNR1 (GPR91) | AequoZen Frozen Cells | ES-744-AF |
| | SUCNR1 (GPR91) | AequoScreen Cell Line | ES-744-A |
| Thyrotropin Releasing Hormone | TRH ₁ | AequoScreen Cell Line | ES-700-A |
| | TRH ₁ | AequoZen Frozen Cells | ES-700-AF |
| Urotensin | (Uro II, GPR 14) -rat | AequoScreen Cell Line | ES-441-A |
| | (Uro II, GPR 14) -rat | AequoZen Frozen Cells | ES-441-AF |
| | UT (Uro II, GPR14) | AequoScreen Cell Line | ES-440-A |
| Vasoactive Intestinal Peptide | PAC ₁ (PACAP) | AequoScreen Cell Line | ES-272-A |
| | VPAC ₁ (VIP) | AequoScreen Cell Line | ES-273-A |
| | VPAC ₂ | AequoScreen Cell Line | ES-274-A |
| Vasopressin | V _{1A} | AequoScreen Cell Line | ES-361-A |
| | V _{1A} | AequoZen Frozen Cells | ES-361-AF |
| | V _{1B} | AequoScreen Cell Line | ES-362-A |
| | V ₂ | AequoScreen Cell Line | ES-363-A |
| Ion Channel Cell Lines | | | |
| 5-Hydroxytryptamine | 5-HT _{3A} | AequoScreen Cell Line | ES-402-A |
| | 5-HT _{3A} | AequoZen Frozen Cells | ES-402-AF |
| Transient receptor potential | TRPA1 | PhotoScreen Cell Line | AX-004-PCL |
| Parental cell lines | | | |
| Parental cell lines | 1321N1 | AequoScreen | ES-000-A27 |
| | 1321N1 Gα ₁₆ | AequoScreen | ES-000-A28 |
| | CHO-K1 | AequoScreen | ES-000-A12 |
| | CHO-K1 Gα ₁₆ | AequoScreen | ES-000-A24 |
| | CHO-K1 Gαq ₁₅ | AequoScreen | ES-000-A21 |
| | HEK293 | AequoScreen | ES-000-A30 |
| | HEK293 Gα ₁₆ | AequoScreen | ES-000-A26 |
| | Parental Frozen Cells | CHO-K1 | AequoZen Frozen Cells |
| CHO-K1 Gα ₁₆ | | AequoZen Frozen Cells | ES-000-A24F |
| Parental cell lines | | | |
| Plasmids | Gα ₁₆ | AequoScreen | ES-003-AC |
| | | AequoScreen | ES-002-AC |

PerkinElmer offers a wide range of luminescence microplate readers and liquid handling solutions to suit your needs. PerkinElmer can easily automate all your luminescence applications. PerkinElmer respects the intellectual property of GPCR patents. IP-protected products may vary from country to country.

THROUGHPUT SOLUTIONS FOR EVERY NEED

PerkinElmer offers a wide range of luminescence microplate readers to suit your needs. From the the new EnSpire with ultra-sensitive luminescence for basic research and assay development, to the MicroBeta² LumiJET platform for drug discovery, PerkinElmer can easily deliver the optimum solution for all your luminescence applications.



EnSpire Multilabel Plate Reader



VICTOR X Light Luminescence
Plate Reader



EnVision Multilabel Plate Reader

EnSpire Multilabel Plate Reader

EnSpire offers affordable detection across multiple platforms including ultra sensitive luminescence, Alpha technology and fluorescence/absorbance using quad-monochromator technology. With ultra-sensitive luminescence, you can use fewer precious cells to obtain the same robust signal as standard luminescence. Extra sensitivity compensates for low transfection assays, enhancing assay viability. The stronger signal increases your assay window, so there is no need to redesign assays. All this combined can save up to 75% in reduced transfection reagents and substrates.

EnSpire leverages EnVision's world class performance technology to deliver flexibility and sensitivity, with easy to use software, operational temperature control and bottom reading capability.

VICTOR™X Light Plate Reader

VICTOR X Light is a flexible microplate-based benchtop reader that offers considerably more sensitivity for luminescence measurement than standard readers—using either aequorin or luciferase assays. It is ideal for a variety of applications including cell-based assays, immunoassays, toxicity screens, and gene expression. Use it as a convenient stand-alone system or, with optional accessories, integrate it into robotic systems as needed. It is ideal for PerkinElmer's mix-and-measure luminescence assay technologies.

VICTOR X Multilabel Reader

With over 10 counting modes, multiple options, and a range of detection technologies, including luminescence, the VICTOR X has the right configuration for your assay needs. It can operate as a stand-alone instrument or integrate into robotic systems. This is your solution if you perform fast multi-wavelength assays in a range of microplates.

EnVision Multilabel Reader

For fast, high performance luminescence detection, consider the EnVision Multilabel Plate Reader. Offered with standard luminescence and optional ultra-sensitive luminescence, this flexible instrument includes two detectors enabling simultaneous dual wavelength reading, below emission reading, barcode readers, a high speed light source and adjustment of measurement height function. Besides luminescence, the EnVision reads all non-isotopic detection technologies, including fluorescence intensity, fluorescence polarization, time-resolved fluorescence (TRF) and absorbance, in all major applications. It's designed to be expanded as your throughput and application requirements increase. The instrument is fully field upgradable with stackers, temperature control, dispenser, and/or AlphaScreen capabilities. EnVision instruments are easily integrated into robotic systems, and accept microplates from 1- to 3456-wells.

MicroBeta² and TopCount Plate Counters

The MicroBeta² and TopCount are the most popular scintillation and luminescence counters in the world for low and medium throughput applications.

Thanks to their high speed photon counting circuitry, MicroBeta², MicroBeta² LumiJET and TopCount have best-in-class glow-type luminescence performance with linear responses in excess of 20 million CPS and backgrounds below 100 CPS. Temperature control features guarantee the optimum counting conditions for luminescence assays.

If your applications require flash luminescence, choose MicroBeta² LumiJET with a choice of single or dual channel injectors. Run dual aequorin (agonist and antagonist) screens with no crosstalk between wells for improved data accuracy and assay performance and no need of an external dispensing device. And with 1, 2, 6 or 12-detector configurations, multiple detectors can be used simultaneously to maximize your throughput in formats up to 384-well plates.

ViewLux Screening Platform

ViewLux is an ultra-high throughput microplate imager for high sensitivity and fast measurement of light from fluorescence polarization (FP), fluorescence intensity, time-resolved fluorescence (TRF), luminescence and absorbance assays.


Because the instrument reads entire plates in one exposure, throughput is not affected by plate density. Using 1536-well plates and, for example, fluorescence intensity, ViewLux allows throughputs in excess of 200,000 samples per hour under continuous operation.

The instrument supports both robot loading and batch mode operation. Up to 64 plates can be loaded for unattended operation. Users can alternate between batch and robot loading according to need.

The JANUS Automated Workstation

From compound management and sample preparation to downstream applications for drug discovery, molecular biology and clinical research, we provide innovative solutions for your automated liquid handling needs. The JANUS Automated Workstation delivers precision liquid handling and flexible pipetting capabilities for pre-dilutions and serial dilutions.

Need a more specialized automation solution? Our Integration Solutions Team stands ready to work with you to develop a custom solution that fits your needs, and our decades of experience mean we have the technical and applications expertise to completely automate your application – maximizing your efficiency and results quality.



MicroBeta² LumiJET
Plate Counter



ViewLux uHF
Microplate Imager



JANUS Automated Workstation

LUMINESCENCE INSTRUMENT QUICK REFERENCE GUIDE

| | | EnSpire | VICTOR X Light |
|------------------------------------|---|--|--|
| FEATURES | Throughput Glow | Up to 50,000 samples per day | Up to 50,000 samples per day |
| | Throughput Flash (Aequorin and Photina) | | Up to 1,000 samples per day (96-well) |
| | Additional/Optional read modes | Dedicated Ultra-sensitive luminescence, FI, Absorbance, Alpha technology | Dedicated Luminescence |
| | Temperature control | +2 °C up to 65 °C, ±0.5 °C | +2 °C to 45 °C |
| | Automation | Robot interface | ActiveX/COM interfaces |
| | Injectors | NA | 1-4 |
| | Injection volumes | NA | 5–350 µL in 1 µL increments |
| | Shaker | Linear, Orbital & Double Orbital | Linear, Orbital & Double Orbital |
| | Microplate format | 1-384 well plates | 1-384 well plates |
| LUMINESCENCE APPLICATIONS | Flash | | X |
| | Glow | X | X |
| | Dual | X | X |
| | Dual emission | | X |
| | Aequorin | | X |
| | Kinetics | X | X |
| LUMINESCENCE SPECIFICATIONS | Bottom Read | X | NA |
| | Detector | Low noise photo multiplier tube | Low noise photo multiplier tube |
| | Dynamic Range | 6 logs | 6 logs |
| | Sensitivity | ATP < 10 pM (Glow) | ATP typically < 10 pM in glow assays |
| | | Firefly luciferase < 1 fg | ATP typically < 100 fM in flash assays |
| | Crosstalk | <0.02% for grey plates | < 0.02% for 96-well B&W Isoplates; < 0.02% for 384-well |

DETECTION INSTRUMENT

| EnVision | ViewLux | MicroBeta ² | MicroBeta ² LumiJET |
|---|-------------------------------------|---------------------------------------|---|
| Up to 100,000 samples per day | 200,000 samples per day | Up to 100,000 samples per day | Up to 100,000 samples per day |
| Up to 1,200 samples per day (96-well) | N/A | | Up to 9,000 samples per day (384-well), 12 detector, 30s per sample |
| Ultra-sensitive Luminescence option, TRF, FI, FP, Absorbance, Alpha technology | TRF, FI, FP, Absorbance | Radiometric Luminescence | Radiometric Luminescence |
| +2 °C to 45 °C | +0 °C to 37 °C | PMT cooling option; no temp on sample | PMT cooling option; no temp on sample |
| Robot loading | Robot loading | Robot loading | Robot loading |
| 1 - 2 | NA | NA | 1 - 2 per detector, up to 12 detectors |
| 2-475 µL in 1 µL increments | NA | NA | 5-250 µL in 1 µL increments |
| Linear, Orbital & Double Orbital | NA | NA | NA |
| 1-3456-well plate | 1-3456-well plates | 24- and 96 OR 96- and 384-well plates | 24- and 96 OR 96- and 384-well plates |
| X | | | X |
| X | X | X | X |
| X | | X | X |
| X | | | |
| X | | | X |
| X | X | X | X |
| X | X | X | X |
| Dedicated low noise PMT for Luminescence, 1-2 parallel PMT's for other technologies | -100 °C cooled CCD | 1, 2, 6 or 12 parallel PMTs | 1, 2, 6 or 12 parallel PMTs |
| 6.5 logs | 3.5 logs | 5.5 logs | 5.5 logs |
| ATP (384-well plate, 50 µL) < 10 pM | ATP (384-well plate, 50 µL) < 10 pM | ATP <10 pM (glow) | ATP <10pM (glow) |
| | | Firefly luciferase <1 fg | Firefly luciferase <1 fg |
| < 0.02% (382-well) | | 0.002% in white 96-well OptiPlate | 0.002% in white 96-well OptiPlate |



APPLICATION-FOCUSED MICROPLATES

OptiPlate™: white polystyrene OptiPlate microplates provide excellent light reflection and the highest efficiency with low background for luminescence. Available in 96-well, 384-well and 1536-well formats.

ProxiPlate™: a shallow well design brings the reagent into closer proximity to the reader's detectors and increases signal. Available in 96-well, 384-well and 1536-well formats.

CulturPlate™: ideal for work with cell-based applications, providing a sterile, tissue culture treated environment with an opaque well bottom. Available in 96-well, 384-well and 1536-well formats.

ViewPlate®: optimized for work with cell-based applications providing a sterile, tissue culture treated environment with a clear well bottom for microscopic visualization or bottom reading. Available in 96-well, 384-well and 1536-well formats.



PerkinElmer microplates for every need.

| Microplate | Color (Frame/Bottom) | Applications | Format (Wells) | Catalog No. | Packaging (case of) |
|-----------------------|----------------------|---|----------------|---------------------------------------|---------------------|
| OptiPlate | white/white | Opaque plates for general luminescence assays and suspension cell-based assays | 96 | 6005291 / 6005290 / 6005299 | 5 / 50 / 200 |
| | | | 384 | 6007291 / 6007290 / 6005299 | 5 / 50 / 200 |
| | | | 1536 | 6004291 / 6004290 / 6004299 | 5 / 50 / 200 |
| ProxiPlate | white/white | Opaque, shallow-well for assay miniaturization | 96 | 6006291 / 6006290 / 6006299 | 5 / 50 / 200 |
| | | | 384 | 6008281 / 6008280 / 6008289 | 5 / 50 / 200 |
| CulturPlate | white/white | Opaque, tissue-culture treated, sterile for cell-based assays and luminescence applications | 96 | 6005681 / 6005680 / 6005688 / 6005689 | 3 / 50 / 160 / 200* |
| | | | 384 | 6007681 / 6007680 / 6007688 / 6007689 | 3 / 50 / 160 / 200* |
| | | | 1536 | 6004681 / 6004680 / 6004688 / 6004691 | 3 / 50 / 160 / 200* |
| ViewPlate | white/clear | Clear bottomed, tissue culture treated, sterile for cell-based luminescence applications | 96 | 6005263 / 6005181 | 3 / 50 |
| | | | 384 | 6007481 / 6007480 | 3 / 40 |
| | | | 1536 | 6004481 / 6004480 | 2 / 40 |
| 1/2 Area Plate | white/white | Opaque plates for general assay miniaturization | 96 | 6005561 / 6005560 / 6005568 | 5 / 50 / 200 |

*Without lids. Order lids separately.

For a complete product listing, please visit:
www.perkinelmer.com/microplates

ACCESS TO
SCIENTISTS BEHIND
THE TECHNOLOGY



OnPointSM services

Time-challenged?

Would added expertise help you advance?

When your research demands outstrip your valuable internal resources, look to PerkinElmer to provide customized solutions to meet your luminescence detection needs.

Our services are tailored, whether you need to outsource the entire assay development process, or simply get some consultative support—we can provide the level of service you need:

- Assay development
- Custom labeling services
- Custom microplate barcoding and coating
- Custom radiosynthesis
- Application/new product development
- Automation and liquid handling solutions
- System integration
- OEM partnerships

Custom GPCR cell line generation

PerkinElmer can generate stably or transiently transfected custom GPCR cell lines for your assay readout of choice, including:

- Competition and saturation ligand binding
- GTP γ S binding assays
- Functional assay readout for cAMP, Inositol Phosphate, ERK phosphorylation Ca²⁺ flux using luminescence photoproteins, fluorescence or both

Ready-to-use frozen cells

We can culture and freeze cells from any of our catalog frozen cells and cell line products using an optimized protocol. They can either be supplied propagatable or irradiated to stop all replication.

- Batch QC criteria (both on EC₅₀ and S/N ratio)
- Validated cells that are ready to use for cAMP or luminescence calcium testing
- Available in off-the-shelf aliquots of 10 million cells+

GPCR assay development

Our scientists can also help you with your customized binding or cellular GPCR assay using your choice of readout.

For a complete description of OnPoint services please visit:

www.perkinelmer.com/onpoint



**For science. For all the possibilities.
For a better tomorrow.**

At PerkinElmer, we share our customers' commitment to finding answers to the mysteries of human health. With proven expertise in reagents, assays, cellular imaging, detection systems, automation and micropates, we offer the right combination of technologies and service to enable scientists around the world to rapidly discover new therapies.

Contact us at 800.762.4000 or (+1) 203.925.4602
or [visit www.perkinelmer.com/contactus](http://www.perkinelmer.com/contactus).

PerkinElmer, Inc.
940 Winter Street
Waltham, MA 02451 USA
P: (800) 762-4000 or
(+1) 203-925-4602
www.perkinelmer.com



For a complete listing of our global offices, visit www.perkinelmer.com/ContactUs

Copyright © 2009-2010, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. All other trademarks are the property of their respective owners.

007091G_01

Printed in USA

Jun. 2010